


EndoC- β H5

Unlimited quantity of ready-to-use functional human beta cells



Bruno BLANCHI, PhD
b.blanchi@humancelldesign.com

Hamza OLLEIK, PharmD, PhD
h.olleik@humancelldesign.com

Glucose dose response similar to native beta cells

Glucose stimulated insulin secretion (GSIS) response :


- Dose dependent

Highest potentiation between **5.5 and 8 mM** Glucose

- Highly sensitive

> **10 fold** increase in insulin secretion

- ✓ **A real physiological human beta cell model**





Reproducible dose dependent responses to GLP-1 and GIP receptor agonists

Exendin-4 and [D-Ala 2]-GIP dose responses :

- Reproducible responses
- Dose dependent
- Robust (> 3.5) fold increase
- Reproducible EC50

✓ Robust pharmacological responses




GSIS reproducibility between batches :

- Reproducible response to glucose
- Reproducible response to GLP-1 receptor agonist (Exendin-4)
- High absolute values of insulin secretion

✓ Batch to batch reproducibility

✓ Systematic QC validation



Glucose (mM) - 20 20
Exendin-4 (nM) - - 1

Take Home Messages: Added Values



EndoC-βH5

- ✓ **Batch to Batch Reproducibility**
- ✓ **Flexibility: Thaw cells when needed**
- ✓ **Availability: Unlimited**
- ✓ **HTS: 96 and 384 well plate compatibility**
- ✓ **Chronic Treatment: >4 weeks**
- ✓ **Time saving: Thaw-and-go - Results in few days**

EndoC-βH5[®]

Independant Validation

6



Article

Large-Scale Functional Genomics Screen to Identify Modulators of Human β-Cell Insulin Secretion

Iwona Szczerbinska ^{1,*}, Annamaria Tessitore ², Lena Kristina Hansson ³, Asmita Agrawal ², Alejandro Ragel Lopez ², Marianne Helenius ^{3,4}, Andrzej R. Malinowski ¹, Barak Gilboa ², Maxwell A. Ruby ¹, Ramneek Gupta ³ and Carina Ämmälä ¹

¹ Department of Discovery Biology and Pharmacology, Novo Nordisk Research Centre Oxford,

Using the EndoC-βH5[®] model, a new publication was published validating novel genes regulating human β-cell insulin secretion.


1. Established a **384-well** functional screen
2. Used **disease-relevant** insulin secretion **endpoints**
3. Successfully applied **siRNA**-mediated loss-of-function
4. Identified potential **T2D drug targets**



ENDOC-BH CELLS
Bruno BLANCHI, PhD
Hamza OLLEIK, PharmD, PhD
All rights reserved © Human Cell Design

Cell Comparison To Native β Cells

		ENDOC- β H1	ENDOC- β H5	NATIVE β CELLS
Functionality	Physiological glucose dose response	+	+++	+++
	GLP-1/GIP response	No	Yes	Yes
	Insulin content (μ g/M C)	0.5 – 1	Up to 10	Up to 10
Characteristics	Proliferation	Yes	No	No
	Fully mature	No	Yes	Yes
	Amplification	> 100 passage	Single use	Single use
	Purity	100% β cells	100% β cells	α / β / δ cells
Experimental	Time for doing your first GSIS	8 weeks	7 days	Islet preparation
	Chronic Treatment	Yes	> 4 weeks	Few days
	Handling	Culture and Preparation	Thaw-and-go / Ready to use	Preparation
	Reproducibility and Robustness	+	+++	+
	Flexibility	+	+++	-
Logistics	Availability	Unlimited	Unlimited	Rare
Screening	96 well plate	No	Yes	Hard due to availability
	HTS	Conditional	Yes	Hard due to availability



HUMAN CELL DESIGN
ACCELERATING DISCOVERY

ENDOC-BH CELLS
Bruno BLANCHI, PhD
Hamza OLLEIK, PharmD, PhD
All rights reserved © Human Cell Design

Acknowledgement



Bruno BLANCHI, PhD
b.blanchi@humancelldesign.com

Hamza OLLEIK, PharmD, PhD
h.olleik@humancelldesign.com



Hamza OLLEIK, PharmD, PhD

h.olleik@humancelldesign.com

www.humancelldesign.com

